

PORTABLE HEAVY TRACE METALS LABORATORIES

HM1000 METALYSER® PORTABLE

An economy level trace heavy metals analyser, offering a total of 10 different parameters down to 3ppb, with completely automated analysis and interpretation on the powerful handheld unit.

Available parameters

As (III), As (III + V), Cd, Cr (VI) Cu, Hg, Mn, Ni, Pb, Zn



Metalometer HM2000



HM2000 METALYSER® DELUXE

A multi-technology instrument, offering electrochemical and photometric techniques to enable testing of a wider range of 13 parameters at a wide spectrum of detection limits.

Available parameters

Al, As (III), As (III + V), B, Cd, Cr (VI) Cu, Fe, Hg, Mn, Ni, Pb, Zn

HM3000 METALYSER® FIELD PRO

An advanced trace heavy metals analyser, with integrated waterproof tablet PC and powerful bespoke software, both of which enable the user to analyse our widest range of metals, at detection limits as low as 0.1ppb.

Available parameters

Au, As (III), As (III + V), Bi, Cd, Co, Cr (VI) Cu, Hg, Mn, Ni, Pb, Sb, Sn, Tl, Zn



